

MINIMED™ 780G SYSTEM

DETAILED MENU MAP

History & Graph

- History
 - Summary
 - Daily History
 - Alarm History
- Sensor Glucose Review
 - High Limit
 - Low Limit
 - Days to Average
- Graph
 - Time in Range

Status

- Suspend All Delivery
- SmartGuard Checklist
- Pump
- Sensor

Insulin

- Bolus
 - Bolus Wizard
 - Manual Bolus
 - Delivery Settings
- Basal
 - Temp Basal
 - Basal Patterns
 - Basal Pattern Setup
- Suspend all Delivery/
Resume Basal
- Delivery settings

Reservoir & Set

- New Reservoir
- Fill Cannula

Paired Devices

- Pair New Device
- Pair CareLink
Mobile
Meter
Sensor

SmartGuard

- SmartGuard Checklist
- Temp Target
- SmartGuard Settings
 - Target
 - Auto Correction
- SmartGuard On and Off

Sound & Vibration

- Silence Sensor Alerts
- Volume
- Sound
- Vibration
- Alert Settings

Blood Glucose

- BG

Settings

- Alert Settings**
 - High Alert
 - Low Alert
 - Snooze High & Low
 - Reminders
- Delivery settings**
 - Bolus Wizard Setup
 - Basal Pattern Setup
 - Max Basal/Bolus
 - Dual/Square Wave
 - Bolus Increment
 - Bolus Speed
 - Preset Bolus Setup
 - Preset Temp Setup
- Device Settings**
 - Sensor On and Off
 - Time & Date
 - Device Info
 - Display
 - Block Mode
 - Self Test
 - Review Settings
 - Manage Settings
 - Easy Bolus
 - Auto Suspend
 - Language

